

BT-200 Vascular Doppler HI-dop

HI-dop is showing you how the blood is flowing.

HI-dop is pocket size doppler for physicians, surgeons and other clinical specialists.

- Compact and light
- High-Sensitivity Doppler probe
- Ergonomic design
- 4 Ranges of probes
- Easy to use
- Long time continuous use (6 hours)



▶ 2 · 4 · 5 · 8MHz Available

BT-200 Vascular Doppler **Hi-dop**

Specifications

Display	<ul style="list-style-type: none"> • LCD Back light • Singnal quality display 	<ul style="list-style-type: none"> • Low battery display
Heart rate range	<ul style="list-style-type: none"> • 50~240 bpm 	
Accuracy	<ul style="list-style-type: none"> • $\pm 2\%$ 	
Ultrasound Frequency	<ul style="list-style-type: none"> • 2, 4, 5, 8 MHz 	
Ultrasound Intensity	<ul style="list-style-type: none"> • $< 10 \text{ m W/cm}^2$ 	
Speaker output power	<ul style="list-style-type: none"> • 1.2 W (Built in Loudspeaker) 	
Audio Output for Earphone		
Auto Shut off Sound	<ul style="list-style-type: none"> • 1 Minute 	
Auto Shut off Power	<ul style="list-style-type: none"> • 5 Minutes 	
Power	<ul style="list-style-type: none"> • 1.5 V Battery X 2 (Type : AA) 	
PC Interface	<ul style="list-style-type: none"> • Sound card 	
Battery life	<ul style="list-style-type: none"> • 360 Minutes 	
Accessories	<ul style="list-style-type: none"> • Carrying case • Gel (60ml) 	<ul style="list-style-type: none"> • Battery (1.5V AA X 2)

