

WONJIN ~ DVT 7700 Compressible Limb Therapy System

DVT-7700 helps stimulate the venous system and address the factors associated with clot formation. This therapy helps increase the volume and velocity of flow from superficial veins up through the deep veins.

- ✓ Assists in decreasing venous stasis, blood pooling, and venous congestion.
- ✓ Helps decrease blood vessel distension and helps lower the risk of injury to the endothelium tissue.
- ✓ Enhances fibrinolysis, which assists in activating anti-clotting factors.
- ✓ Improve blood flow in the patient's legs, reducing the risk for blood clots.



Technical Specifications:

| | |
|--------------------------|---|
| Power requirement | 100V~240VAC, 50/60Hz |
| Power Consumption | Maximum 18W |
| Power Cord | Hospital Grade Plug (Detachable) |
| Operation Mode | A1,A2,A3,B1,B2,B3,B4,C1,D1,D2 |
| Adjustable Timer | 0-99 Hours |
| Applied Pressure | Leg sleeve : 30~60 mmHg (Calf, Thigh) Foot Cuffs : 130mmHg |
| Operation Mode | Continuous |
| Compression Type | Leg sleeve: Gradient, Sequential |
| Dimension | 200mm (L) * 190mm (W) * 260 mm (H) |
| Weight | 2.3 Kg |
| Battery | Capacity 10.8 Vdc 4,400 mAh Lithiumion (optional) Run Time 6-8 Hours |

Sold and Serviced by:

BIO VISION MEDICAL SYSTEMS

No:49&50, Second floor, Samuvel Nagar
Vadaperumbakkam, Chennai – 600 060.
Ph No: 73388 88672/044-4868 3042
E-mail: sales@biovisionchennai.com
Web: www.biovisionchennai.com

