

# endo E-CUBE

NEW



**S**mart design  
**A**dvanced  
**E**asy

## Features

- Auto Reverse function (E CUBE-Black)
- Torque control from 0.6 to 5.2 N.cm
- gear ratio : 16 : 1
- 9 Programmable preset memorizing
- Foot control On/Off switch
- Up to 5-hour battery life
- Compact designed controller including various function
- Supply with stable torque
- operation available with battery or adapter

## Product Specification

### Control Unit

Rated input	DC20V 0.5A
Output	DC4.5V 0.4A
Charging time	5 hours approx.
Dimensions	W 110 X D 196 X H 139mm
Weight	582g

### Motor (CUBE EM)(CUBE EP)


Rated input	DC4.5V 0.4A
RPM	13,000rpm

### Angle Handpiece (ACL(B)-42EP)

Rated input	16 : 1
Dimensions	φ 16.7 X L63mm
Weight	34.8g

## Product Components

- (1) Control Unit (2type)
- (2) AC Adapter
- (3) Contra Angle ( 1EA )
- (4) Motor ( 1EA )
- (5) Handpiece Stand
- (6) On/Off Switch (Option)



- A type On/Off Button
- Contra Angle Option ----- ● B type Non On/Off Button

