


Introducing brilliance to the oral cavity!!

LUMINRG

Illuminated Mirror – Two in One



- Bright, focused light
- Different LED colors for multiple functions
- Autoclave Safe
- Compatible with most standard dental mirrors and instruments
- Adjustable angle and distance from LED
- Long lasting LED
- Easily replaceable standard AAA battery
- Ergonomic and easy to use
- Attractive pricing



**SOLD BY / UNDER AMC WITH
BIO VISION MEDICAL SYSTEMS**

Old # 46, (New No. 105), 1st Floor, St. Mary's Road,
Abhiramapuram, Chennai - 600 018. Tel : 2461 8393

MEDICNRG™
Precise Endo Technology


Optimizing Your Endo Performance

Introducing brilliance to the oral cavity!!

LUMINRG

Illuminated Mirror – Two in One

Introduces a new level of bright light inside the oral cavity and provides dental professionals with a whole new level of comfort.



Brilliant Light: The simplicity of the device is brilliant, literally and metaphorically. The bright LED light introduces a new level of light into the dark oral cavity, giving the dental professionals a whole new level of comfort. A bright, focused and miniature light which can reduce the need to use the bulky overhead lights.

Multiple Functions: The LumiNRG can perform different functions, by using the different colored LEDs that can be interchanged. In order to diagnose craze lines, fractures and cracks in the tooth, attach the blue LED. The bright white LED or the orange LED can support curing processes. The white LED for faster curing, the orange LED for slower curing.

Adaptable: LumiNRG is designed to be used with standard threaded or unthreaded mirrors. Based on a new patent which enables the adjustment of the mirror to any angle or depth from the light.


Versatile: In addition to use with standard mirrors, the LumiNRG can also be used with many different tips and instruments, providing they are of a standard 3 mm threaded diameter. Scalers, pluggers and other tips can easily be attached to the LumiNRG. The desired instrument's angle can be adjusted for optimizing the light source.

Autoclavable! The LumiNRG's illumination head can be placed directly into the autoclave along with the other dental tools. There is no need to detach the mirror or other instrument prior to sterilization. The battery pack must be detached.

Long Life: The battery consumption is regulated by a smart electronic circuit, enabling the battery to last for hundreds of treatments with a constant light intensity. Uses only 1 standard AAA battery.

Ergonomic: Say goodbye to the bulky and inconvenient overhead light. No longer will dental professionals have to adjust and re-adjust the large overhead light during the procedure.

Minimal Size: The LumiNRG ensures that virtually no extra space will be taken up in the already limited space of the oral cavity and is especially helpful when working on the remote molar teeth.



MedicNRG Europe

Tel: +31-165-30 86 86,
E-mail: info@medicnrg.com/nl
Web: www.medicnrg.com/eu

MedicNRG USA

Toll-free: +1 (888) 429-0240,
E-mail: infoUSA@medicnrg.com
Web: www.medicnrg.com/usa

Analytical & Medical Technologies

Sun Palm View, Plot No 1-1A,
Unit No B-703, Sector 15,
Sanpada, Navi Mumbai - 400705
Maharashtra, INDIA.
Tel.: 022-6514 1277 / 093226 73077
Email: amtcorp@sify.com
Website: www.amtcorp.in

MEDICNRGTM
Precise Endo Technology