

MediJECT - The Professional's Choice

**INTRA LIGAMENTARY SYRINGE**



For

- a. Local Anesthesia of a Single upper or lower tooth
- b. Pedodontia-avoids numbness and therefore a possible trauma to lip, tongue or cheek
- c. Treatment in several quadrants possible-important for instance when extracting a bicuspid on each side for orthodontic reasons in same setting
- d. Treatment of a single tooth-for instance an occlusal obturation, especially at the beginning of a caries in a bicuspid or molar
- e. Infected root canals and other endodontic treatments.
- f. Deep curettage and periodontal- surgical interventions e Professional's Choice



**Details:**

- Unique ratchet design for smooth delivery of anesthesia
- Removable headpiece with angled threaded nozzle
- Servo-closing lever for optimal control and support for index finger
- Multiplies your strength nine times
- 

**Features:**

- Automatic limitation of the power pressure
- No back pressure
- Single tooth anesthesia
- Supports regional block anesthesia
- Eliminates risk of haematoma
- Autoclavable

**BIO VISION MEDICAL SYSTEMS**  
**Old # 45 (New #103), 1<sup>st</sup> Floor, St. Mary's Road, Abhiramapuram,**  
**Chennai – 600 018& Tel: 2461 8393 Telefax: 91-44-24618393**  
**E-mail: [biovisionmedicalsistemas@gmail.com](mailto:biovisionmedicalsistemas@gmail.com)**