

Diamond DeluxTM

because we care

An ISO 9001 : 2008 & 13485 : 2012 and CE Certified Company

COMFORT 60

MULTI PURPOSE
ELECTRIC OPERATING TABLES
* WITH LONGITUDINAL TOP SLIDE



ISO
9001

ISO
13485

CE

■ Technical Specifications

Table top size (LxW)	: 1930 x 530 mm
Min height (without mattress)	: 730 mm
Max height (without mattress)	: 980 mm
Trendelenburg / Reverse Trendelenburg	: 20°/20°
Lateral Tilt (Left/Right)	: 20°/20°
Back section adjustment	: +70°/-25°
Head section adjustment	: +15°/-90°
Leg section adjustment	: 90° DOWN & SPLIT
Load capacity	: 180 Kg.
Power supply	: 220 VAC -50Hz

■ Standard Accessories

'L' Shaped Anaesthetic Screen	: 1 Pc.
Padded Shoulder Support	: 1 Pair
Padded Lateral Support	: 1 Pair
Radio-translucent Armboard with pad	: 2 Pcs.
Goepel Knee Crutch with straps	: 1 Pair
PU Foam Mattress	: 1 Set

